



YO Proteins AB Karolinska Institutet Science Park
Alfred Nobels allé 10, 141 52 Huddinge, Sweden

www.yoproteins.com

Heparin – agarose resin

Description

Heparin (pig) immobilized on 4% agarose (slurry)

Catalog Number 346

Amount 1 ml

Format

50% slurry (2.0 ml), gel (1.0 ml) in 0.2 M NaCl, 0.1 M Na-borate, 0.005 M Na-azide, pH 8.0

Ligand

Heparin (pig)

Labeling

5 mg heparin per 1 ml of resin

Matrix

4% cross-linked agarose

Binding Capacity

6 mg of protein per 1 ml of packed gel

Storage

Shipping at ambient temperature. Storage at +4°C for 12 months. Do not freeze

RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, DRUG OR
HOUSEHOLD USE